

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 372 053

Micro burette, 5 ml, str. stopc., PTFE, AS, conf. cert.




Batch: 19-22

Serial:

Measuring Nr.: 58710

Batch size: 9

Volume:	5,0000 ml	3,0000 ml	1,0000 ml
Tolerance:	± 0,0100 ml	± 0,0100 ml	± 0,0100 ml
Average value:	4,9956 ml	2,9943 ml	0,9994 ml
Deviation:	0,0016 ml	0,0006 ml	0,0009 ml
Precision (CV):	0,0311 %	0,0204 %	0,0905 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	08.09.2022
Balance:	Sartorius - MSE224S-0CE -DU

DE-M