

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 372 053

Micro burette, 5 ml, str. stopc., PTFE, AS, conf. cert.




Batch: 12-23

Serial:

Measuring Nr. 69700

Batch size: 6

Volume:	5,0000 ml	3,0000 ml	1,0000 ml
Tolerance:	± 0,0100 ml	± 0,0100 ml	± 0,0100 ml
Average value:	4,9977 ml	2,9957 ml	0,9972 ml
Deviation:	0,0037 ml	0,0024 ml	0,0005 ml
Precision (CV):	0,0745 %	0,0809 %	0,0528 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	77
Date:	23.10.2024
Balance:	Sartorius - LE244S-OCE

DE-M