

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 372 103

Micro burette, 10 ml, str.stopc., PTFE, AS, conf. cert.




Batch: 43-20

Serial:

Measuring Nr. 43326

Batch size: 6

Volume:	10,000 ml	6,000 ml	2,000 ml
Tolerance:	± 0,0200 ml	± 0,0200 ml	± 0,0200 ml
Average value:	10,0014 ml	5,9932 ml	1,9954 ml
Deviation:	0,0072 ml	0,0060 ml	0,0069 ml
Precision (CV):	0,0723 %	0,0995 %	0,3437 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	9
Date:	15.12.2020
Balance:	Sartorius - MSE224S-0CE -DU

DE-M