

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 383 022

Micro bur., Bang, 2 ml, lat. stopcock, AS, conf. cert.

Batch: 09-16




Serial:

Measuring Nr. 12105

Batch size:

1

Volume:	2,0000 ml	1,2000 ml	0,4000 ml
Tolerance:	± 0,0100 ml	± 0,0100 ml	± 0,0100 ml
Measured value:	2,0039 ml	1,1986 ml	0,3999 ml



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

4

Date:

18.07.2017

Balance:

Sartorius -
LE244S-OCE

DE-M