

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 383 022

**Micro bur., Bang, 2 ml, lat. stopcock, AS, conf. cert.**




**Batch:** 32-18

**Serial:**

**Measuring Nr.** 34953

**Batch size:** 2

| Volume:                | 2,0000 ml   | 1,2000 ml   | 0,4000 ml   |
|------------------------|-------------|-------------|-------------|
| <b>Tolerance:</b>      | ± 0,0100 ml | ± 0,0100 ml | ± 0,0100 ml |
| <b>Average value:</b>  | 2,0004 ml   | 1,2021 ml   | 0,3986 ml   |
| <b>Deviation:</b>      | 0,0025 ml   | 0,0035 ml   | 0,0011 ml   |
| <b>Precision (CV):</b> | 0,1242 %    | 0,2897 %    | 0,2729 %    |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

|                  |                           |
|------------------|---------------------------|
| <b>Operator:</b> | 9                         |
| <b>Date:</b>     | 09.12.2019                |
| <b>Balance:</b>  | Sartorius -<br>LE244S-OCE |

