

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 383 025

**Micro burettes, 2 ml, lat. needle valve, AS, conf. cert.**

**Batch:** 32-18




**Serial:**

**Measuring Nr.** 38354

**Batch size:**

8

Volume:	2,0000 ml	1,2000 ml	0,4000 ml
Tolerance:	± 0,0100 ml	± 0,0100 ml	± 0,0100 ml
Average value:	1,9995 ml	1,2023 ml	0,3983 ml
Deviation:	0,0041 ml	0,0023 ml	0,0031 ml
Precision (CV):	0,2039 %	0,1901 %	0,7803 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

**Operator:**

9

**Date:**

02.06.2020

**Balance:**

Sartorius -  
LE244S-OCE

DE-M