

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 383 025

Micro burettes, 2 ml, lat. needle valve, AS, conf. cert.




Batch: 41-20

Serial:

Measuring Nr.: 56508

Batch size: 10

Volume:	2,0000 ml	1,2000 ml	0,4000 ml
Tolerance:	± 0,0100 ml	± 0,0100 ml	± 0,0100 ml
Average value:	2,0010 ml	1,2000 ml	0,3966 ml
Deviation:	0,0009 ml	0,0023 ml	0,0019 ml
Precision (CV):	0,0472 %	0,1890 %	0,4695 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	07.06.2022
Balance:	Sartorius - MSE224S-0CE -DU

