

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 383 052

Micro burette Bang, 5 ml, lat. stopc., AS, conf.cert.

Batch: 40-18




Serial:

Measuring Nr. 27518

Batch size:

2

Volume:	5,0000 ml	3,0000 ml	1,0000 ml
Tolerance:	± 0,0100 ml	± 0,0100 ml	± 0,0100 ml
Average value:	5,0063 ml	3,0027 ml	1,0021 ml
Deviation:	0,0017 ml	0,0030 ml	0,0000 ml
Precision (CV):	0,0344 %	0,0988 %	0,0018 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	9
Date:	06.02.2019
Balance:	Sartorius - LE244S-OCE

DE-M