

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 383 055

Micro burettes, 5 ml, lat. needle valve, AS, conf. cert.




Batch: 40-18

Serial:

Measuring Nr. 25360

Batch size: 10

Volume:	5,0000 ml	3,0000 ml	1,0000 ml
Tolerance:	± 0,0100 ml	± 0,0100 ml	± 0,0100 ml
Average value:	5,0023 ml	2,9998 ml	0,9981 ml
Deviation:	0,0035 ml	0,0034 ml	0,0027 ml
Precision (CV):	0,0697 %	0,1120 %	0,2698 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	05.11.2018
Balance:	Sartorius - LE244S-OCE

