

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 515 910

Graduated pipettes class AS 25 ml




Batch: 16-19

Serial:

Measuring Nr.: 29278

Batch size: 500

Volume:	25,000 ml	15,000 ml	5,000 ml
Tolerance:	± 0,1000 ml	± 0,1000 ml	± 0,1000 ml
Average value:	25,0005 ml	15,0033 ml	5,0142 ml
Deviation:	0,0324 ml	0,0234 ml	0,0134 ml
Precision (CV):	0,1298 %	0,1562 %	0,2675 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	16.04.2019
Balance:	Sartorius - LE244S-OCE

