

# Certificate of performance

**Article:** 4711632

**Burettes, 50 ml, lateral valve, Schellb., AS, conf. cert.**




**Batch:** 43-20

**Serial:**

**Measuring Nr.:** 43401

**Batch size:** 50

Volume:	50,000 ml	30,000 ml	10,000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	49,9982 ml	30,0054 ml	10,0287 ml
Deviation:	0,0158 ml	0,0163 ml	0,0070 ml
Precision (CV):	0,0317 %	0,0542 %	0,0698 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	9
Date:	21.12.2020
Balance:	Sartorius - MSE224S-0CE -DU