

Certificate of performance



ISO 9001
CERTIFIED

Article: 4711632

Burettes, 50 ml, lateral valve, Schellb., AS, conf. cert.




Batch: 17-22

Serial:

Measuring Nr. 58458

Batch size: 28

Volume:	50,000 ml	30,000 ml	10,000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	50,0113 ml	29,9992 ml	10,0072 ml
Deviation:	0,0151 ml	0,0217 ml	0,0109 ml
Precision (CV):	0,0302 %	0,0724 %	0,1090 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	31.08.2022
Balance:	Sartorius - LE244S-OCE

