

Certificate of performance



ISO 9001
CERTIFIED

Article: 4711632

Burettes, 50 ml, lateral valve, Schellb., AS, conf. cert.




Batch: 31-22

Serial:

Measuring Nr. 61812

Batch size: 50

Volume:	50,000 ml	30,000 ml	10,000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	50,0155 ml	30,0279 ml	10,0300 ml
Deviation:	0,0054 ml	0,0088 ml	0,0046 ml
Precision (CV):	0,0108 %	0,0295 %	0,0462 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

7

Date:

01.03.2023

Balance:

Sartorius -
LE244S-OCE

DE-M