

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 4711632

**Burettes, 50 ml, lateral valve, Schellb., AS, conf. cert.**




**Batch:** 11-23

**Serial:**

**Measuring Nr.** 69497

**Batch size:** 27

| Volume:         | 50,0000 ml  | 30,0000 ml  | 10,0000 ml  |
|-----------------|-------------|-------------|-------------|
| Tolerance:      | ± 0,0500 ml | ± 0,0500 ml | ± 0,0500 ml |
| Average value:  | 50,0033 ml  | 30,0098 ml  | 10,0039 ml  |
| Deviation:      | 0,0120 ml   | 0,0181 ml   | 0,0124 ml   |
| Precision (CV): | 0,0239 %    | 0,0602 %    | 0,1239 %    |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

|           |                                   |
|-----------|-----------------------------------|
| Operator: | 4                                 |
| Date:     | 15.10.2024                        |
| Balance:  | Sartorius -<br>MSE224S-0CE<br>-DU |

