

Certificate of performance



ISO 9001
CERTIFIED

Article: 4711632

Burettes, 50 ml, lateral valve, Schellb., AS, conf. cert.




Batch: 01-26

Serial:

Measuring Nr. 77303

Batch size: 60

Volume:	50,0000 ml	30,0000 ml	10,0000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	50,0126 ml	30,0032 ml	9,9993 ml
Deviation:	0,0107 ml	0,0148 ml	0,0153 ml
Precision (CV):	0,0213 %	0,0493 %	0,1527 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	09.04.2026
Balance:	Sartorius - MSE224S-0CE -DU

