

Certificate of performance



ISO 9001
CERTIFIED

Article:

4711632

Burettes, 50 ml, lateral valve, Schellb., AS, conf. cert.

Batch:

18-26

Serial:

Measuring Nr.

78134

Batch size:

30

Volume:

50,0000 ml

30,0000 ml

10,0000 ml

Tolerance:

$\pm 0,0500$ ml

$\pm 0,0500$ ml

$\pm 0,0500$ ml

Average value:

49,9826 ml

29,9925 ml

9,9919 ml

Deviation:

0,0093 ml

0,0073 ml




0,0098 ml

Precision (CV):

0,0186 %

0,0244 %

0,0985 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

7

Date:

27.05.2026

Balance:

Sartorius -
LE244S-OCE

DE-M