

Certificate of performance



ISO 9001
CERTIFIED

Article: 4712333

Pellet-burettes, 50ml, lat. needle valve, AS, conf. cert.

Batch: 30-17




Serial:

Measuring Nr. 22078

Batch size:

12

Volume:	50,000 ml	30,000 ml	10,000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	49,9990 ml	30,0254 ml	10,0216 ml
Deviation:	0,0089 ml	0,0032 ml	0,0141 ml
Precision (CV):	0,0178 %	0,0108 %	0,1407 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

4

Date:

09.07.2018

Balance:

Sartorius -
LE244S-OCE

DE-M