

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 715 000

**Graduated cylinder, 5:0,1 ml, A, conf. cert.**




**Batch:** 48-24

**Serial:**

**Measuring Nr.** 70846

**Batch size:** 800

Volume:	1,0000 ml	3,0000 ml	5,0000 ml
<b>Tolerance:</b>	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
<b>Average value:</b>	0,9884 ml	2,9982 ml	4,9991 ml
<b>Deviation:</b>	0,0139 ml	0,0122 ml	0,0101 ml
<b>Precision (CV):</b>	1,3915 %	0,4053 %	0,2024 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	4
<b>Date:</b>	02.12.2024
<b>Balance:</b>	Sartorius - MSE224S-0CE -DU

