

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 715 001 BL USP

**Grad. cylinder, 10:0,2 ml, blue grad., USP, A, conf. cert.**




**Batch:** 04-19

**Serial:**

**Measuring Nr.:** 26979

**Batch size:** 100

Volume:	2,0000 ml	6,0000 ml	10,0000 ml
Tolerance:	± 0,1000 ml	± 0,1000 ml	± 0,1000 ml
Average value:	1,9972 ml	6,0024 ml	10,0023 ml
Deviation:	0,0233 ml	0,0169 ml	0,0142 ml
Precision (CV):	1,1636 %	0,2810 %	0,1416 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	9
Date:	24.01.2019
Balance:	Sartorius - LE244S-OCE

