

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 001 BL

Grad. cylinders, 10:0,2 ml, blue grad, A, conf. cert.




Batch: 04-19

Serial:

Measuring Nr. 27030

Batch size: 1200

Volume:	2,0000 ml	6,0000 ml	10,0000 ml
Tolerance:	± 0,1000 ml	± 0,1000 ml	± 0,1000 ml
Average value:	1,9925 ml	6,0069 ml	10,0241 ml
Deviation:	0,0202 ml	0,0178 ml	0,0208 ml
Precision (CV):	1,0121 %	0,2959 %	0,2076 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	28.01.2019
Balance:	Sartorius - LE244S-OCE

