

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 001

Grad. cylinder, 10:0,2 ml, A, conf.cert.




Batch: 38-17

Serial:

Measuring Nr. 14681

Batch size: 500

Volume:	2,0000 ml	6,0000 ml	10,0000 ml
Tolerance:	± 0,1000 ml	± 0,1000 ml	± 0,1000 ml
Average value:	2,0011 ml	6,0092 ml	10,0173 ml
Deviation:	0,0184 ml	0,0130 ml	0,0106 ml
Precision (CV):	0,9178 %	0,2169 %	0,1063 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	05.10.2017
Balance:	Sartorius - LE244S-OCE

