

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 001

Grad. cylinder, 10:0,2 ml, A, conf.cert.




Batch: 32-20

Serial:

Measuring Nr. 39695

Batch size: 280

Volume:	2,0000 ml	6,0000 ml	10,0000 ml
Tolerance:	± 0,1000 ml	± 0,1000 ml	± 0,1000 ml
Average value:	2,0083 ml	6,0049 ml	10,0208 ml
Deviation:	0,0201 ml	0,0182 ml	0,0179 ml
Precision (CV):	1,0044 %	0,3035 %	0,1790 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	26.08.2020
Balance:	Sartorius - LE244S-OCE

