

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 001

Grad. cylinder, 10:0,2 ml, A, conf.cert.




Batch: 14-26

Serial:

Measuring Nr. 77556

Batch size: 499

| Volume: | 2,0000 ml | 6,0000 ml | 10,0000 ml |
|-----------------|-------------|-------------|-------------|
| Tolerance: | ± 0,1000 ml | ± 0,1000 ml | ± 0,1000 ml |
| Average value: | 1,9828 ml | 6,0125 ml | 10,0211 ml |
| Deviation: | 0,0114 ml | 0,0093 ml | 0,0115 ml |
| Precision (CV): | 0,5703 % | 0,1556 % | 0,1149 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

| | |
|-----------|---------------------------|
| Operator: | 7 |
| Date: | 20.04.2026 |
| Balance: | Sartorius - LE244S-OCE |

