

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 715 003 BL USP

**Grad. cyl., 25:0,5 ml, blue grad., USP, A, conf. cert.**




**Batch:** 10-19

**Serial:**

**Measuring Nr.** 28245

**Batch size:** 50

Volume:	5,0000 ml	15,0000 ml	25,0000 ml
Tolerance:	± 0,1700 ml	± 0,1700 ml	± 0,1700 ml
Average value:	5,0330 ml	15,0039 ml	24,9779 ml
Deviation:	0,0314 ml	0,0372 ml	0,0284 ml
Precision (CV):	0,6275 %	0,2477 %	0,1135 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	4
Date:	11.03.2019
Balance:	Sartorius - LE244S-OCE

DE-M