


Certificate of performance

Article number	3 715 003 BL	Graduated Cylinder 25 ml Class AS KB			
Batch	03-17	Volume	5,0000ml	15,0000ml	25,0000ml
Serial number		Tolerance	0,2500ml	0,2500ml	0,2500ml
Measuringnr./Batch size	10060 / 279	Average volume	4,9965 ml	15,0039 ml	25,0223 ml
Balance	Sartorius - LE244S-0CE	Deviation	0,0420 ml	0,0261 ml	0,0189 ml
Operator	4	Number of measurings	15	15	15
Date	20.01.2017				



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
 97877 Wertheim
 Germany
 +49 (0) 9342 / 9301-0
 | info@witeg.de

