

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 715 003 BL

**Grad. cylinder, 25:0,5 ml, blue grad., A, conf. cert.**




**Batch:** 03-18

**Serial:**

**Measuring Nr.** 17757

**Batch size:** 280

Volume:	5,0000 ml	15,0000 ml	25,0000 ml
<b>Tolerance:</b>	± 0,2500 ml	± 0,2500 ml	± 0,2500 ml
<b>Average value:</b>	4,9939 ml	15,0187 ml	25,0361 ml
<b>Deviation:</b>	0,0524 ml	0,0445 ml	0,0453 ml
<b>Precision (CV):</b>	1,0478 %	0,2968 %	0,1811 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	9
<b>Date:</b>	01.02.2018
<b>Balance:</b>	Sartorius - LE244S-OCE

