

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 003 BL

Grad. cylinder, 25:0,5 ml, blue grad., A, conf. cert.




Batch: 34-23

Serial:

Measuring Nr. 64154

Batch size: 280

Volume:	5,0000 ml	15,0000 ml	25,0000 ml
Tolerance:	± 0,2500 ml	± 0,2500 ml	± 0,2500 ml
Average value:	4,9736 ml	14,9850 ml	25,0096 ml
Deviation:	0,0435 ml	0,0369 ml	0,0240 ml
Precision (CV):	0,8708 %	0,2463 %	0,0958 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	26.08.2023
Balance:	Sartorius - MSE224S-0CE -DU

DE-M