

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 003 BL

Grad. cylinder, 25:0,5 ml, blue grad., A, conf. cert.




Batch: 43-25

Serial:

Measuring Nr. 75011

Batch size: 150

Volume:	5,0000 ml	15,0000 ml	25,0000 ml
Tolerance:	± 0,2500 ml	± 0,2500 ml	± 0,2500 ml
Average value:	5,0247 ml	15,0298 ml	25,0320 ml
Deviation:	0,0494 ml	0,0420 ml	0,0545 ml
Precision (CV):	0,9881 %	0,2798 %	0,2181 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	27.10.2025
Balance:	Sartorius - LE244S-OCE

