

Certificate of performance



Article: 3 715 003 USP

Graduated cylinder, 25:0,5 ml, USP, A, conf. cert.




Batch: 43-25

Serial:

Measuring Nr. 75033

Batch size: 150

Volume:	5,0000 ml	15,0000 ml	25,0000 ml
Tolerance:	± 0,1700 ml	± 0,1700 ml	± 0,1700 ml
Average value:	5,0405 ml	15,0389 ml	25,0425 ml
Deviation:	0,0372 ml	0,0305 ml	0,0277 ml
Precision (CV):	0,7435 %	0,2036 %	0,1107 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	27.10.2025
Balance:	Sartorius - LE244S-OCE

