

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 715 003

**Graduated cylinder, 25:0,5 ml, A, conf. cert**




**Batch:** 43-25

**Serial:**

**Measuring Nr.** 75068

**Batch size:** 600

Volume:	5,0000 ml	15,0000 ml	25,0000 ml
<b>Tolerance:</b>	± 0,2500 ml	± 0,2500 ml	± 0,2500 ml
<b>Average value:</b>	4,9740 ml	14,9966 ml	25,0168 ml
<b>Deviation:</b>	0,0383 ml	0,0313 ml	0,0261 ml
<b>Precision (CV):</b>	0,7656 %	0,2087 %	0,1044 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	7
<b>Date:</b>	28.10.2025
<b>Balance:</b>	Sartorius - LE244S-OCE

