

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 005 BL USP

Grad. cylinder, 50:1 ml, blue grad., USP, A, conf. cert.




Batch: 28-21

Serial:

Measuring Nr.: 47890

Batch size: 100

Volume:	10,000 ml	30,000 ml	50,000 ml
Tolerance:	± 0,2500 ml	± 0,2500 ml	± 0,2500 ml
Average value:	9,9850 ml	30,0722 ml	50,1170 ml
Deviation:	0,0621 ml	0,0519 ml	0,0589 ml
Precision (CV):	0,6207 %	0,1729 %	0,1178 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	16.07.2021
Balance:	Sartorius - LE244S-OCE

