

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 005 BL

Graduated cylinders, 50 ml, blue grad., A, conf. cert.




Batch: 44-17

Serial:

Measuring Nr. 15700

Batch size: 280

Volume:	10,000 ml	30,000 ml	50,000 ml
Tolerance:	± 0,5000 ml	± 0,5000 ml	± 0,5000 ml
Average value:	10,0492 ml	30,0565 ml	50,0758 ml
Deviation:	0,0544 ml	0,0363 ml	0,0348 ml
Precision (CV):	0,5438 %	0,1209 %	0,0696 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	14.11.2017
Balance:	Sartorius - LE244S-OCE

