

# Certificate of performance



**Article:** 3 715 005 BL

**Graduated cylinders, 50:1 ml, blue grad., A, conf. cert.**




**Batch:** 05-21

**Serial:**

**Measuring Nr.** 44375

**Batch size:** 500

Volume:	10,000 ml	30,000 ml	50,000 ml
Tolerance:	± 0,5000 ml	± 0,5000 ml	± 0,5000 ml
Average value:	10,0395 ml	30,0253 ml	50,0368 ml
Deviation:	0,0689 ml	0,0511 ml	0,0528 ml
Precision (CV):	0,6887 %	0,1704 %	0,1057 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	7
Date:	15.02.2021
Balance:	Sartorius - LE244S-OCE

