

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 005 BL

Graduated cylinders, 50:1 ml, blue grad., A, conf. cert.




Batch: 46-23

Serial:

Measuring Nr.: 65238

Batch size: 499

Volume:	10,000 ml	30,000 ml	50,000 ml
Tolerance:	± 0,5000 ml	± 0,5000 ml	± 0,5000 ml
Average value:	9,9338 ml	29,9773 ml	49,9963 ml
Deviation:	0,0581 ml	0,0375 ml	0,0254 ml
Precision (CV):	0,5807 %	0,1250 %	0,0507 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	22.11.2023
Balance:	Sartorius - MSE224S-0CE -DU

