

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 715 005 USP

**Grad.cylinder, 50:1 ml, USP, A, conf.cert.**




**Batch:** 40-20

**Serial:**

**Measuring Nr.:** 41229

**Batch size:** 100

Volume:	10,000 ml	30,000 ml	50,000 ml
Tolerance:	± 0,2500 ml	± 0,2500 ml	± 0,2500 ml
Average value:	9,9163 ml	29,9822 ml	50,0276 ml
Deviation:	0,0479 ml	0,0284 ml	0,0260 ml
Precision (CV):	0,4789 %	0,0947 %	0,0520 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

**Operator:**

4

**Date:**

**Balance:**

Sartorius -  
MSE224S-0CE  
-DU

