

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 005 USP

Grad.cylinder, 50:1 ml, USP, A, conf.cert.




Batch: 28-21

Serial:

Measuring Nr. 47889

Batch size: 280

Volume:	10,0000 ml	30,0000 ml	50,0000 ml
Tolerance:	± 0,2500 ml	± 0,2500 ml	± 0,2500 ml
Average value:	9,9677 ml	30,0337 ml	50,0980 ml
Deviation:	0,0854 ml	0,0747 ml	0,0465 ml
Precision (CV):	0,8541 %	0,2491 %	0,0929 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	15.07.2021
Balance:	Sartorius - LE244S-OCE

