

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 715 005 USP

**Grad.cylinder, 50:1 ml, USP, A, conf.cert.**




**Batch:** 51-21

**Serial:**

**Measuring Nr.** 51300

**Batch size:** 280

Volume:	10,0000 ml	30,0000 ml	50,0000 ml
Tolerance:	± 0,2500 ml	± 0,2500 ml	± 0,2500 ml
Average value:	10,0371 ml	30,0257 ml	50,0527 ml
Deviation:	0,0487 ml	0,0517 ml	0,0306 ml
Precision (CV):	0,4868 %	0,1725 %	0,0613 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	7
Date:	12.01.2022
Balance:	Sartorius - LE244S-OCE

