

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 005 USP

Graduated cylinder, 50:1 ml, USP, A, conf.cert.




Batch: 47-22

Serial:

Measuring Nr. 60750

Batch size: 150

Volume:	10,000 ml	30,000 ml	50,000 ml
Tolerance:	± 0,2500 ml	± 0,2500 ml	± 0,2500 ml
Average value:	9,9633 ml	30,0073 ml	50,0770 ml
Deviation:	0,0499 ml	0,0449 ml	0,0441 ml
Precision (CV):	0,4989 %	0,1496 %	0,0881 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	16.12.2022
Balance:	Sartorius - MSE224S-0CE -DU

