

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 005 USP

Graduated cylinder, 50:1 ml, USP, A, conf.cert.




Batch: 17-24

Serial:

Measuring Nr. 66969

Batch size: 150

| Volume: | 10,000 ml | 30,000 ml | 50,000 ml |
|------------------------|-------------|-------------|-------------|
| Tolerance: | ± 0,2500 ml | ± 0,2500 ml | ± 0,2500 ml |
| Average value: | 9,9890 ml | 29,9644 ml | 49,8948 ml |
| Deviation: | 0,0507 ml | 0,0239 ml | 0,0184 ml |
| Precision (CV): | 0,5069 % | 0,0795 % | 0,0368 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

| | |
|------------------|---------------------------|
| Operator: | 7 |
| Date: | 25.04.2024 |
| Balance: | Sartorius - LE244S-OCE |

