

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 715 005 USP

**Graduated cylinder, 50:1 ml, USP, A, conf.cert.**




**Batch:** 50-24

**Serial:**

**Measuring Nr.** 71222

**Batch size:** 280

Volume:	10,0000 ml	30,0000 ml	50,0000 ml
<b>Tolerance:</b>	± 0,2500 ml	± 0,2500 ml	± 0,2500 ml
<b>Average value:</b>	9,9751 ml	29,9842 ml	49,9738 ml
<b>Deviation:</b>	0,0708 ml	0,0754 ml	0,0689 ml
<b>Precision (CV):</b>	0,7080 %	0,2514 %	0,1378 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	4
<b>Date:</b>	13.12.2024
<b>Balance:</b>	Sartorius - MSE224S-0CE -DU

