

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 005 USP

Graduated cylinder, 50:1 ml, USP, A, conf. cert.




Batch: 47-25

Serial:

Measuring Nr. 75942

Batch size: 280

Volume:	10,0000 ml	30,0000 ml	50,0000 ml
Tolerance:	± 0,2500 ml	± 0,2500 ml	± 0,2500 ml
Average value:	10,0188 ml	30,0628 ml	50,0397 ml
Deviation:	0,0366 ml	0,0286 ml	0,0348 ml
Precision (CV):	0,3664 %	0,0954 %	0,0697 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	10.12.2025
Balance:	Sartorius - MSE224S-0CE -DU

