

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 005 USP

Graduated cylinder, 50:1 ml, USP, A, conf. cert.




Batch: 12-26

Serial:

Measuring Nr. 77235

Batch size: 149

Volume:	10,0000 ml	30,0000 ml	50,0000 ml
Tolerance:	± 0,2500 ml	± 0,2500 ml	± 0,2500 ml
Average value:	10,0024 ml	29,9810 ml	49,8864 ml
Deviation:	0,0739 ml	0,0727 ml	0,0465 ml
Precision (CV):	0,7391 %	0,2425 %	0,0930 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator: 7

Date: 01.04.2026

Balance: Sartorius -
LE244S-OCE

DE-M