

Certificate of performance



Article: 3 715 005

Graduated cylinders, 50:1 ml, A, conf. cert.




Batch: 23-23

Serial:

Measuring Nr.: 63374

Batch size: 1183

Volume:	10,000 ml	30,000 ml	50,000 ml
Tolerance:	± 0,5000 ml	± 0,5000 ml	± 0,5000 ml
Average value:	9,9904 ml	30,0262 ml	50,0172 ml
Deviation:	0,0759 ml	0,0468 ml	0,0440 ml
Precision (CV):	0,7590 %	0,1560 %	0,0880 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	16.06.2023
Balance:	Sartorius - LE244S-OCE

