

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 005

Graduated cylinders, 50:1 ml, A, conf. cert.




Batch: 46-23

Serial:

Measuring Nr.: 65199

Batch size: 1200

Volume:	10,0000 ml	30,0000 ml	50,0000 ml
Tolerance:	± 0,5000 ml	± 0,5000 ml	± 0,5000 ml
Average value:	9,9636 ml	29,9997 ml	50,0219 ml
Deviation:	0,0657 ml	0,0485 ml	0,0436 ml
Precision (CV):	0,6565 %	0,1615 %	0,0872 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	20.11.2023
Balance:	Sartorius - LE244S-OCE

