

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 010 BL

Grad. cylinder, 100:1,0 ml, blue grad., A, conf. cert.




Batch: 18-22

Serial:

Measuring Nr. 55321

Batch size: 280

Volume:	20,0000 ml	60,0000 ml	100,0000 ml
Tolerance:	± 0,5000 ml	± 0,5000 ml	± 0,5000 ml
Average value:	19,8967 ml	59,9605 ml	100,0610 ml
Deviation:	0,0874 ml	0,0892 ml	0,1202 ml
Precision (CV):	0,4372 %	0,1486 %	0,1202 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	10.05.2022
Balance:	Sartorius - LE244S-OCE

