

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 010 USP

Graduated cylinders, 100:1 ml, USP, A, conf. cert.




Batch: 19-26

Serial:

Measuring Nr. 78062

Batch size: 500

Volume:	20,0000 ml	60,0000 ml	100,0000 ml
Tolerance:	$\pm 0,5000$ ml	$\pm 0,5000$ ml	$\pm 0,5000$ ml
Average value:	20,0962 ml	60,0997 ml	100,0628 ml
Deviation:	0,0591 ml	0,0671 ml	0,0611 ml
Precision (CV):	0,2955 %	0,1119 %	0,0611 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

7

Date:

13.05.2026

Balance:

Sartorius -
LE244S-OCE

DE-M