

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 715 010

**Graduated cylinder, 100:1,0 ml, A, conf. cert.**




**Batch:** 38-22

**Serial:**

**Measuring Nr.** 59314

**Batch size:** 2095

Volume:	20,000 ml	60,000 ml	100,000 ml
<b>Tolerance:</b>	± 0,5000 ml	± 0,5000 ml	± 0,5000 ml
<b>Average value:</b>	20,1090 ml	60,1248 ml	100,1922 ml
<b>Deviation:</b>	0,1644 ml	0,1099 ml	0,1102 ml
<b>Precision (CV):</b>	0,8218 %	0,1831 %	0,1102 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	4
<b>Date:</b>	22.09.2022
<b>Balance:</b>	Sartorius - LA 1200 S

