

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 715 010

**Graduated cylinder, 100:1,0 ml, A, conf. cert.**




**Batch:** 18-26

**Serial:**

**Measuring Nr.** 78245

**Batch size:** 1527

|                        |                   |                   |                    |
|------------------------|-------------------|-------------------|--------------------|
| <b>Volume:</b>         | <b>20,0000 ml</b> | <b>60,0000 ml</b> | <b>100,0000 ml</b> |
| <b>Tolerance:</b>      | $\pm 0,5000$ ml   | $\pm 0,5000$ ml   | $\pm 0,5000$ ml    |
| <b>Average value:</b>  | 20,0825 ml        | 60,0811 ml        | 100,0711 ml        |
| <b>Deviation:</b>      | 0,0926 ml         | 0,0578 ml         | 0,0634 ml          |
| <b>Precision (CV):</b> | 0,4628 %          | 0,0963 %          | 0,0634 %           |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

|                  |                                   |
|------------------|-----------------------------------|
| <b>Operator:</b> | 4                                 |
| <b>Date:</b>     | 10.06.2026                        |
| <b>Balance:</b>  | Sartorius -<br>MSE224S-0CE<br>-DU |

DE-M