

# Certificate of performance



**Article:** 3 715 025 BL

**Grad. cylinder, 250:2 ml, blue grad., A, conf. cert.**




**Batch:** 50-18

**Serial:** -

**Measuring Nr.:** 26558

**Batch size:** 280

Volume:	50,0000 ml	150,0000 ml	250,0000 ml
<b>Tolerance:</b>	± 1,0000 ml	± 1,0000 ml	± 1,0000 ml
<b>Average value:</b>	50,0243 ml	149,9294 ml	249,7398 ml
<b>Deviation:</b>	0,1891 ml	0,1240 ml	0,1130 ml
<b>Precision (CV):</b>	0,3782 %	0,0827 %	0,0452 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	4
<b>Date:</b>	18.12.2018
<b>Balance:</b>	Sartorius - LA 1200 S

