

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 025 USP

Graduated cylinders, 250:2 ml, A, USP, conf. cert.




Batch: 44-23

Serial:

Measuring Nr. 64912

Batch size: 150

Volume:	50,000 ml	150,000 ml	250,000 ml
Tolerance:	± 1,000 ml	± 1,000 ml	± 1,000 ml
Average value:	50,0291 ml	150,1092 ml	250,1384 ml
Deviation:	0,1894 ml	0,2282 ml	0,2520 ml
Precision (CV):	0,3788 %	0,1521 %	0,1008 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	06.11.2023
Balance:	Sartorius - LE244S-OCE

